



ENDOSCOPE SERIES SURGICAL INSTRUMENTS



LAPAROSKOPİ CERRAHİ EL ALETLERİ



Yenişehir Mah. Turgut Özal Bulvarı No: 36 D: 14 Ataşehir / İstanbul

+90 216 455 80 40

info@yetkinmedikal.com www.yetkinmedikal.com

EndoGlobe

Laparoscopic Instruments

Laparoscope

0° Laparoscope
Φ10×330mm
Φ5×330mm



30° Laparoscope
Φ10×330mm
Φ5×330mm



Trocar and Sealing Cap

flip type

801.008.1 Φ12.5 × 95mm
801.004.1 Φ10.5 × 95mm
801.003.1 Φ5.5 × 95mm



flip type(with autoshield safety)

801.163.1 Φ12.5 × 95mm
801.075.1 Φ10.5 × 95mm
801.076.1 Φ5.5 × 95mm



cross-shaped

801.008 Φ12.5 × 95mm
801.004 Φ10.5 × 95mm
801.003 Φ5.5 × 95mm
804.002 Φ3.5 × 60mm



cross-shaped(with autoshield safety)

801.163 Φ12.5 × 95mm
801.075 Φ10.5 × 95mm
801.076 Φ5.5 × 95mm



trocar pyramidal

trocar with
autoshield,safety

magnetic type

801.008.3 Φ12.5 × 95mm
801.004.3 Φ10.5 × 95mm
801.003.3 Φ5.5 × 95mm



magnetic type(with autoshield safety)

801.163.3 Φ12.5 × 95mm
801.075.3 Φ10.5 × 95mm
801.076.3 Φ5.5 × 95mm



flip type

801.003.2 Φ5 × 95mm
801.004.2 Φ10 × 95mm



flip type(with autoshield safety)

801.076.2 Φ5 × 95mm
801.076.2 Φ5 × 95mm



trocar cannula
without metal thread

trocar cannula
with metal thread

trocar cambered

802.008,trocar with
metal thread,Φ5.5mm



802.009,trocar with
metal thread,Φ8.5mm



802.010,trocar with
metal thread,Φ10.5mm



802.011,trocar with
metal thread,Φ12.5mm



801.108,new threaded
trocar Φ5.5×88mm




801.109,new threaded
trocar Φ10.5×88mm



801.086.2, Φ10.5mm sealing cap/washer



801.086.1, Φ5.5mm sealing cap/washer



801.086, Φ3mm sealing cap/washer

Reducer

801.010 sleeve type, Φ10.5-Φ5.5mm



801.010 100mm
801.010.1 120mm
801.010.2 150mm
801.010.3 170mm



801.010.4 22mm



801.145 bayonet type
Φ10.5-Φ5.5mm



801.146 screw type
Φ10.5-Φ5.5mm

Aspiration needle

801.027 Φ5 × 330mm



801.027.1 7#

801.027.4 18#

801.027.2 9#

801.027.5 20#

801.027.3 16#

801.027.6 22#

Bile duck Knife

801.025

801.025.1

801.026

801.025 Φ5 × 330mm

801.025.1 Φ5 × 330mm

801.026 Φ10 × 330mm

Monopolar HF Electrode

Φ5 × 330mm



808.017 monopolar
hook electrode "L"

808.017.3 monopolar
hook electrode "J"

808.019.1 monopolar
electric spade

808.017.1 monopolar
sharp "L" hook

808.021.1 monopolar
electric knife

808.020.1 monopolar
electric pin

808.017.2 monopolar
electric hook

808.018.1 monopolar
coagulation

Veress Needle

801.011 $\Phi 2.2 \times 110\text{mm}$
801.011.1 $\Phi 2.2 \times 120\text{mm}$
801.011.2 $\Phi 2.2 \times 150\text{mm}$



801.231 Veress Needle Connector

Suction & Irrigation

$\Phi 5/\Phi 10 \times 330\text{mm}$

801.079 push type
801.079.1 push type (big valve)
801.079.2 push type (small valve)



$\Phi 5/\Phi 10 \times 330\text{mm}$

801.149 spring type



5mm Golden Scissors

801.033.1JB
Golden Scissors

$\Phi 5 \times 330\text{mm}$



Bipolar

$\Phi 5 \times 330\text{mm}$

808.044 Bipolar Forceps (screwed)



808.044 doubleaction (straight)



808.044.1 doubleaction (curved)

$\Phi 5 \times 330\text{mm}$

808.045 Bipolar Forceps (three partition)



808.045 single action (straight)

808.045.1 double action (straight)



808.045.2 single action (curved)


808.045.3 double action (curved)



808.045.2 single action (blunt shaped)



808.045.5 single action (round shaped)



Cord for Monopolar



808.082 Cord for Monopolar $\Phi 4 \times 3000\text{mm}$

808.082.1 match OLYMPUS $\Phi 8\text{mm}$

808.082.2 ordinary connect or $\Phi 6\text{mm}$

808.082.6 match HUTONG $\Phi 6.3\text{mm}$

Cord for Bipolar



808.039 Cord for Bipolar 3000mm (match three partition handle)




808.039.1 Cord for Bipolar (match blue handle)


Clip Applier, Polymer Ligation Clip



801.113, XL, $\Phi 10 \times 330\text{mm}$, matches 805.621



801.113.1,L, $\Phi 10 \times 330\text{mm}$, matches 805.622



801.113.2 ,S, $\Phi 5.5 \times 330\text{mm}$, matches 805.623



805.621, XL



805.622, L



805.623, S



804.114.1, Articulated, $\Phi 10 \times 330\text{mm}$, L, matches 805.622

804.114, Articulated, $\Phi 10 \times 330\text{mm}$, L, matches 805.621

Open Titanium Clip Applier




Open Polymer Clip Applier



Clip Remover




Titanium Clip, Clip Applier



Purse Stitching Clamp



Artery/Vein Clip, Clip Applier




Knot Pusher




Clip Applier



Needle Holder



Closure Forceps




Three Leaves Retractor and Five Leaves Retractor



Distraction Forceps(three leaves retractor)



Spring Grasping Forcep



Biliary Tract-Lithotomy Forcep



Golden Finger



Vein/Nerve Retractor




Live Retractor



Bendable forceps



Cholangiography forceps



Φ5×330mm Forceps and Scissors

Handle Φ5mm



801F.001 insulated handle (two partition)



801F.002 insulated handle with ratchet (two partition)



801F.003 insulated handle (three partition)



801F.004 insulated handle (three partition, with ratchet)



801F.013 insulated handle



801F.011 insulated handle

	801.032 straight scissors (single action)		801.133 claw forceps
	801.032.1 straight scissors (double action)		801.037 Maryland forceps 17mm
	801.033 curved scissors (single action)		801.037.1 Maryland forceps 12mm
	801.033.1 curved scissors (double action)		801.037.2 Maryland forceps 14mm
	801.033.2 micro-curved scissors (double action)		801.037.3 Maryland forceps 18mm
	801.034 hook scissors		801.037.4 Maryland forceps 20mm
	801.035 upwards scissors		801.037.5 Maryland forceps 25mm
	801.033.1B insulated scissors		801.037.6 blunt tip Maryland forceps 16mm
	801.036 straight grasper		801.037.7 Maryland forceps (45°)
	801.051 dolphin nose grasper		801.037.8 Maryland forceps (70°)

	801.037.9 Maryland forceps (70°) Mesh tooth		801.195.1 grasper 1+2teeth, double action
	801.038 Mixer right angle		801.031 knot forceps
	801.038.1 Mixer right angle (wider head)		801.044.1 reverse-teeth forceps
	801.038.3 ligature clamp		801.047 cholecystic forceps
	801.037.1 insulated maryland forceps		801.047.1 nine teeth cholecystic forceps
	801.049 fenestrated grasping forceps		801.054 biopsy forceps sharp cupped jaw
	801.048.1 grasping forceps		801.105 claw forceps 2+4teeth, single action
	801.048.2 grasping forceps		801.105.1 claw forceps 2+4teeth, double action
	801.052 appendix grasping forceps		801.044.2 reverse-teeth forceps



801F.005 stainless steel handle



801F.006 stainless steel(self-locking)



801F.008 stainless steel(with upper ratchet)



801F.007 stainless steel handle (with down ratchet)

	801.031.1 B-shaped knot forceps		801.050 bullet nose forceps
	801.062.1 V-shaped grasping forceps		801.057 fallopian tube grasping forceps
	801.061 O-shaped grasping forceps		801.042 duck-jaw grasping forceps
	801.198.1 stomach grasping forceps, curved		801.042.1 duck-jaw grasping forceps coarse teeth
	801.134 babcock grasping forceps		801.043 Maxi grasper
	801.198 stomach forceps (double action)		801.044 reverse-teeth forceps
	801.055 biting&cutting forceps		801.046.1 grasper 2+3 teeth (double action)
	801.057.1 Oviduct forceps		801.058 mucous membrane forceps hollow
	801.053 cupped biopsy forceps		801.058.2 mucous membrane forceps

	801.059 hollow erect grasping forceps		801.030 knot forceps S shape
	801.059.1 hollow erect grasping forceps, long tip		801.061 O-shaped grasping forceps
	801.060.1 Allis grasper longer tip		801.056.1 Debakey grasping forceps(hollow)
	801.062 V-shaped grasping forceps		801.056 Debakey grasping forceps
	801.116 ureteral grasping forceps		801.072.2 sector forceps
	801.045 organ-anchoring grasper		801.056.3 bowel forceps
	801.048 atraumatic forceps		801.056.4 intestinal grasping forceps
	801.049.1 fenestrated grasping forceps		801.018.1 Wolf claws grasping forceps
	801.049.2 endo(lap)clinch fenestrated grasper		

Φ10mm×330mm Forceps and Scissors

Handle(Φ10mm)



801F.010 insulated handle
(two partition)



801F.012 insulated handle
(three partition)

	801.012.2 biopsy forceps matches 801F.010		801.041.1 Mixer dissecting forceps matches 801F.010
	801.039.2 straight scissor matches 801F.010		801.040.2 curved scissors matches 801 F.010
	801.013.2 claw forceps only matches 801F.010		801.012 biopsy forceps matches 801F.016
	801.041 dissecting forceps(right angle) matches 801F.010		801.135 grasping forceps(eight claws) 801F.017
	801.041 dissecting forceps(right angle) matches 801F.010		802.007 (duck-jaw grasping forceps) matches 801F.017
	801.013 claw forceps matches 801F.017		801.014 atraumatic forceps matches 801F.017



801F.015 stainless steel handle, can not be detached by hands(self-locked)




801F.016 stainless steel handle, can not be detached by hands




801F.017 stainless steel handle, can not be detached by hands(with down ratchet)

	801.069 stomach grasping forceps round tip, matches 801F.017		802.002 lobe grasping forceps matches 801F.017
	801.071 anastomotic and tractive grasping forceps matches 801F.017		801.130 US Babcock grasping matches 801F.017
	801.126 reverse teeth matches 801F.017		801.107 gut grasping forceps bendable, matches 801F.017
	802.004 toothless grasping forceps (with wide mouth) matches 801F.017		802.005 grasping forceps (with long chin) matches 801F.017
	801.121 oval grasping forceps matches 801F.017		801.235 interdiction vessel forceps, matches 801F.017
	801.068 stomach grasping forceps matches 801F.017		801.070 suture grasping forceps, matches 801F.017

OC-804 Miniature Instruments for children and Pediatric surgery




Instrument FOR SILS/Single Port Instruments



Detachable Picture of Instruments

Two partition handle



Three partition handle



Assembly Instructions Applicable to Table 1

1. Press button 1, keep handle left part still, and pull the right part 2 till button 1 hidden inside.
2. Slide insert into handle, notice the tongue and groove configuration. The insert should snap firmly into the shaft.
3. Close handle part 2, till you hear a 'click' sound, make sure the insert fits the handle completely.

Disassembly Instructions Applicable to Table 1

1. Press button 1, keep handle left part still, and pull the right part 2 till button 1 hidden inside.
2. Keep pulling 2, till the insert slide out of the handle.
3. Pull out of the insert from the handle.

Assembly Instructions Applicable to Table 2

1. Plug in insert into shaft, keep insert upside and make sure it fits shaft slot perfectly, rotate it clockwise and keep the jaw close.
2. Pull handle 2 part to it's MAX position.
3. Insert shaft into handle till the end of insert touches the handle.
4. Close handle slightly till you hear a 'click' sound.

Disassembly Instructions Applicable to Table 2

1. Press 'PRESS' button, keep handle left part still, and pull the right part 2.
2. Keep pulling 2, till the insert slide out of the handle.
3. Keep jaws open and upside, rotate counterclockwise, and slide insert from the shaft.

Thyroid Surgical Instruments




801.152 separation instruments for subcutaneous tissue $\Phi 10 \times 350\text{mm}$




801.154 bendable detacher $\Phi 5 \times 350\text{mm}$




801.153 thyroid detacher $\Phi 5 \times 350\text{mm}$



801.154 special wire retractor (left curved) $\Phi 1.8/100-250\text{mm}$



801.154 special wire retractor (right curved) $\Phi 1.8/100-250\text{mm}$



801.155 hypodermic injector $\Phi 2.5 \times 250\text{mm}$



801.159 trocar $\Phi 11 \times 125\text{mm}$



801.158 trocar $\Phi 6 \times 125\text{mm}$